



Cigniti

Understanding Disease (COVID) Progression using Chest X-Rays

A Tech Brief



Introduction

COVID-19 is now a global pandemic, with close to 5 million infections globally. As patients get admitted to hospitals, it becomes important for clinicians to understand which of these patients have the potential to turn critical, so that hospitals can plan and prepare accordingly. Longitudinal Chest X-Ray analysis will enable hospitals to prioritize patient care based on the prognostic and severity predictions given by the model. We attempt to provide this insight by generating X-Rays for future time-periods for a patient by studying the previous X-Rays of that patient. This helps the clinician better understand the spread of the disease in a patient, thereby planning better interventions.

In the context of COVID treatment, a deep learning-based CT segmentation of pulmonary opacities is in use. It helps improve the quantification of the disease. Based on U-Net architecture, a convolution neural network is developed to predict expert segmentation.

Our Approach

When there is a multiscale system, we cannot always define the function that governs the system. Whenever possible, it is desirable to work with the linear dynamics of the form.

$$dx/dt = Ax$$

The solution of the above equation is given by:

$$x(t_0 + t) = e^{At} x(t_0)$$

The dynamics are entirely characterized by the eigenvalues and eigenvectors of the matrix A, given by the eigen-decomposition of A:

When A has n distinct eigenvalues, then Λ is a diagonal matrix containing the eigenvalues λ_j and T is a matrix whose columns are the linearly independent eigenvectors ξ_j associated with eigenvalues λ_j . In this case, it is possible to write $A = T\Lambda T^{-1}$, and the solution becomes

$$x(t_0 + t) = T e^{\Lambda t} T^{-1} x(t_0)$$

DMD is a technique to obtain reduced order models for high dimensional systems. It's a purely data-driven method and does not require any knowledge of the underlying equations, originally introduced by Peter Schmidt in fluid dynamics.

As far as we know, this is the first time DMD (Dynamic Mode Decomposition) is being used to study and predict Chest X-Rays in the COVID-19 scenario.

We can extract, from this data, the spatial-temporal coherent structures or patterns that dominate the mentioned data from the dynamical system. DMD is defined to be the eigen decomposition of an approximating best fit linear operator.

Refer to Fig 1 to understand the steps followed in computing standard DMD.

Algorithm 1 (Standard DMD).

1. Arrange the data $\{z_0, \dots, z_m\}$ into matrices

$$X \triangleq [z_0 \ \dots \ z_{m-1}], \quad Y \triangleq [z_1 \ \dots \ z_m]. \quad (2)$$

2. Compute the (reduced) SVD of X (see [32]), writing

$$X = U\Sigma V^*, \quad (3)$$

where U is $n \times r$, Σ is diagonal and $r \times r$, V is $m \times r$, and r is the rank of X.

3. Define the matrix

$$\tilde{A} \triangleq U^* Y V \Sigma^{-1}. \quad (4)$$

4. Compute eigenvalues and eigenvectors of \tilde{A} , writing

$$\tilde{A} w = \lambda w. \quad (5)$$

5. The DMD mode corresponding to the DMD eigenvalue λ is then given by

$$\hat{\varphi} \triangleq U w. \quad (6)$$

6. If desired, the DMD modes can be scaled in a number of ways, as described in Appendix A.

Figure 1: Standard DMD mode reference – [1]

When we feed the data, as it evolves in time, into DMD, we are actually stacking the data into tall vectors evolving in time into vectors; one is shifted delta t in time X and Y. X is any system that exists in high dimensional space n. All the measurements of the vector in n dimensions have n-1 degrees of freedom. There are dominant patterns that are present in this data, so therefore we get spatial-temporal modes and a linear dynamical system.

This means that we get a best fit linear operator A (please refer Fig. 1) that fits and best maps X to Y (please refer Fig. 1). DMD approximates the leading eigen decomposition of the A matrix. The dominant eigenvalues and eigenvectors of A can be reshaped into Eigen flow fields and used to see how this structure evolves in time. It's very much like PCA, where we try to find a low-rank structure and approximate it with a linear dynamical system.

Data

The COVID-19 images and longitudinal information were collected from Dr Cohen’s site [5] using the metadata file.

Results

The data in Figure 2 depicts the X-Rays taken on admission and the following two consecutive days. Using this longitudinal data, the tall X matrix is constructed with a time difference of each passing day.



Figure 2: Original X Rays

After building the Eigen decomposition, each eigenvalue in μ tells us something about the dynamic behaviour of its corresponding DMD mode. If the eigenvalue has a non-zero imaginary part, then there is an oscillation in the corresponding DMD mode. If the eigenvalue is inside the unit circle, then the mode is decaying. If the eigenvalue is outside, then the mode is growing. If the eigenvalue falls exactly on the unit circle, then the mode either grows or decays. From Figure 3, we can observe that one of the modes is growing while the other two are decaying.




Figure 3: Eigen plot

Equipped with the Eigendecomposition of A and a basic understanding of the nature of the system $Y=AX$, it is possible to construct a matrix Ψ . This corresponds to the system’s time evolution in time T . We have used this data to calculate the DMD and were able to predict the data one time-step ahead in time (Figure 4). The fourth X-ray in Fig. 4 was reconstructed from the initial condition and it is our prediction, one step in time for this patient.




Figure 4: Reconstructed Modes

Future Work

We are planning to use the Koopman operator to better understand & estimate the non-linear modes too. We hope to publish this subsequent work in the next couple of weeks.

References

- [1]Kutz, J. Nathan. Data-driven modeling& scientific computation: methods for complex systems & big data. OUP Oxford, 2013.
- [2]Tu, Jonathan H., et al. "On dynamic mode decomposition: theory and applications." arXiv preprint arXiv:1312.0041 (2013). <http://cwwrowley.princeton.edu/papers/Tu-DMD.pdf>
- [3] <https://pubs.rsna.org/doi/full/10.1148/ryct.2020200082>
- [4] https://en.wikipedia.org/wiki/Dynamic_mode_decomposition
- [5] <https://github.com/ieee8023/covid-chestxray-dataset>
- [6] <http://www.pyrunner.com/weblog/2016/07/25/>
- [7] <http://dmdbook.com/>

Success Stories

Senior IT leaders share how our services helped them win in the platform age.

CTO Speak



Data pipeline buildout, Software development, Salesforce development, AWS System admin/DevOps, BI/Dashboard. The execution has been very good.

- Satyadeep "Bobby" Patnaik, CTO



CEO Speak



They understood that product development was iterative and they patiently worked through our requirements even as they rapidly evolved.

- Dr. Ganesh Naidoo, CEO



CTO Speak



Proven technical ability in both web and mobile development; strong project/product management expertise; the ability to become part of the extended BA365 team.

- Graeme Dollar, CTO



COO Speak



I have worked with hundreds of service providers and consultants, RoundSqr (Part of Cigniti) is absolutely one of the best. The people are highly skilled, very hard-working, and have a "can do" attitude.

- Mark Mortimer, COO



Analyst Recognitions



NelsonHall recognized Cigniti as a "LEADER" in its 2022 NEAT evaluation for Overall **Quality** Engineering, Continuous Testing, Application **Security Testing**, and **AI & Cognitive**



Cigniti is mentioned as a "Pure Play Testing Vendor" in Gartner's Market Guide for Application Testing Services, 2022

Cigniti is mentioned as "API Testing Vendor" in Gartner's Hype Cycle for Managed IT Services and APIs, 2022



Recognized as "LEADER" in ISG IPL for Next-Gen ADM Services under Continuous Testing Specialists category for the US region for 2021 and 2022

Recognized as "RISING STAR" in UK Region in ISG IPL for Next-Gen ADM Services 2022



Cigniti is recognized as a **Strong Performer** in the Forrester Wave: **Continuous Automation and Testing** Services, Q3 2021.

One of the top 3 leaders in Agile Testing and DevOps in the Forrester Wave: Continuous Automation and Testing Services, Q3 2017.



Provides Cigniti with "**Best in Class**" rating for Buyer satisfaction.

This tech brief is issued for information only. Cigniti declines all responsibility for any errors and any loss or damage resulting from use of the contents of this tech brief. Cigniti also declines responsibility for any infringement of any third party's Intellectual Property rights but will gladly acknowledge any IPR and correct any infringement of which it is informed.

About Cigniti

Cigniti Technologies Limited (NSE: CIGNITITEC; BSE: 534758) is the World's Leading AI & IP-led Digital Assurance and Digital Engineering Services Company providing software quality engineering, software testing, automation, and consulting services. 4100+ Cignitians worldwide help Fortune 500 & Global 2000 enterprises across 24 countries accelerate their digital transformation journey across various stages of digital adoption and help them achieve market leadership by providing transformation services leveraging IP & Platform-led innovation with expertise across multiple verticals and domains.

Our global customers have benefitted with measurable outcomes, millions of dollars of savings, significant ROI, and delightful, frictionless experiences utilizing our flagship digital assurance and full cycle software quality engineering services. Our AI-led digital engineering services cover Data engineering services, software platform, and digital product engineering, AI/ML engineering services, intelligent automation, big data analytics, and Blockchain development.

We are headquartered in Hyderabad, India, with global offices spread across the USA, Canada, UK, UAE, Australia, South Africa, Czech Republic, and Singapore.

To learn more, visit www.cigniti.com

